

MONTE CARLO STUDY OF THE OPTIMIZATION OF THE MINIBEAM COLLIMATORS IN PROTON MINIBEAM RADIATION THERAPY

C. Guardiola*, C. Peucelle, Y. Prezado

*IMNC – UMR 8165, IN2P3-CNRS, Université Paris 7,
Université Paris 11, Orsay, France*

**guardiola@imnc.in2p3.fr*

INTRODUCTION

AIM: optimizing the minibeam generation in a new radiotherapy approach called **proton-minibeam radiation therapy (pMBRT)**



- Once of the current RT challenges: RT limited by the tolerance doses of normal tissues
- MBRT: Use of a spatial fractionation of the dose & submillimetric field sizes. MBRT has showed increase the tolerance of healthy tissue.
- Pioneer project that offers an optimized way of using protons for therapy.
- Application field: pediatric oncology and treatment of very radioresistant tumors.
- The technique is being technically implemented at Proton Therapy Center Orsay

METHOD: Minibeam Radiation Therapy



Conventional
RT



Spatially fractionated
RT




- Submillimetric field sizes (25 to 700 μm)
- Interbeam separation (400 to 3500 μm)
- Dose profiles: pattern of peaks and valleys



Spatial fractionation

Instead of homogeneous distributions
→ gain in healthy tissue recovery & increase of tolerances [1,2]



Dose-volume effects


The smaller the field size is, the higher the tolerance
→ exponential increase of healthy tissue tolerances [3]

[1] Y. Prezado *et al.*, *J. Synchrotron Radiat.* **19** 60-65 (2012)

[2] P. Deman *et al.*, *Int. J. Radiat. Oncol. Biol. Phys.* **82** 693-700 (2012)


[3] Zeman *et al.*, *Science* (1959)

METHOD: proton MiniBeam Radiation Therapy (pMBRT)




$$PVDR = \frac{D_{\text{peak}}}{D_{\text{valley}}}$$

PVDR has to be as high as possible, while low valley doses are required in order to ensure the preservation of normal tissue



- Gain in healthy tissue recovery & increase of tolerances
- Spatial fraction of the dose in the normal tissue beyond the Bragg peak
- Quasi-homogeneous dose distribution at the Bragg peak location due to multiple Coulomb scattering in depth
- Potential gain in Therapeutic Index




RESULTS: Monte Carlo Simulations



- Monte Carlo simulations (**GATE v7.1**) were used as a method to evaluate the dose distributions of pMBRT irradiations in several configurations of mechanical collimation.
- **Assessed:** the peak and valley depth dose curves, PVDR values, penumbras, and secondary neutron contributions as a function of depth.

Initial configuration:

- 100 MeV proton beams
- 5 multislit collimator (400 μm microslit width, 3200 μm c-t-c)
- 5 cm collimator thickness
- 7 cm phantom-collimator distance



MONTE CARLO SIMULATIONS: Collimator material




- PVDR@W values ~ double other collimators
- Gain in penumbra values: narrower than in conventional radiosurgery (> 2.5 mm) → it helps in normal tissue sparing

MONTE CARLO SIMULATIONS: Neutron contribution

Material	Multislit Collimators(%)	Conventional Collimator (%)
Brass	14.7 ± 0.4	11.5 ± 0.2
Tungsten	39.6 ± 0.9	30.78 ± 0.14
Iron	13.5 ± 0.3	10.6 ± 0.3
Nickel	9.4 ± 0.3	7.4 ± 0.2

- Tungsten delivers the highest secondary neutron yield
- pMBRT delivers a lightly higher secondary neutron yield than conventional collimators ($2 \times 2 \text{ cm}^2$)




D_{neutron} for all depths are
< 0.00016% & 0.0003% of D_{peak}
< 0.0014% & 0.005% of D_{valley}
by using the Brass & W collimators

Biological neutron doses in the phantom (patient) will still be less than 1% of total absorbed dose in the worst case (entrance).

MONTE CARLO SIMULATIONS:


Collimator thickness



if \uparrow collimator thickness: \uparrow gain PVDR, but \downarrow output factor


\rightarrow the number of neutrons produced in the collimator will increase by 25% to deposit the same peak dose

center-to-center distance



if \uparrow c-t-c distance: \uparrow gain PVDR, but \downarrow the homogenization in the target \rightarrow the use of larger c-t-c significantly increases the PVDR, but to obtain a homogeneous dose in the target the interlacing of two orthogonal arrays would be required


MONTE CARLO SIMULATIONS: phantom-collimator distance




If ↓ PCD distance:

- ↓FWHM & ↑PVDR
- ↑Neutron dose, e.g. reducing PCD distance from 7 to 2 cm increases the neutron dose in the surface by a factor of 7.

Full width at half maximum @ Brass



Penumbra @ Brass



CONCLUSIONS

- Tungsten and brass offer the best compromise among the materials evaluated.
- Tungsten multislit provides the highest PVDR and lowest penumbra but it significantly increases the neutron production in the collimator, although the biological doses in the phantom remain below 1%.
- The PVDR are lower for the brass, but it is more advantageous in terms of neutron contamination, manufacturing cost and microetching.
- The aforementioned results will be used to guide the manufacture of optimized collimators for the next pre-clinical studies.

Acknowledgement: Support by the ITMO Cancer AVIESAN 277 (Alliance Nationale pour les Sciences de la Vie et de la Santé, National Alliance for Life Sciences and Health) within the framework of the Cancer Plan (2014-2019). The authors also acknowledge the calculation time granted at Centre de Calcul de Lyon (IN2P3) and at Grand Equipement National de Calcul Intensif (in particular at the supercomputer Curie of CEA) and MareNostrum supercomputational center (Barcelona) is acknowledged.