

**Numerical Calculation of Temperature
Increase During MRI due to the Presence
of Total Hip Replacement Implants**

Vasileios Gkyrinis and Theodoros Samaras
Department of Physics, Aristotle University of Thessaloniki (GR)

Background

- Active or passive implants can heat the surrounding tissue when inside the MRI scanner
- This leads to exclusion of patients from MRI scanning due to implants, unless heating can be assessed and mitigated



Objective

- To demonstrate the possibility of assessing worst-case temperature rise from implants with numerical methods


Temperature increase of tissue due to hip implants inside an MRI birdcage coil
V. Gkyrinis and T. Samaras, Aristotle University of Thessaloniki (GR)

Materials and methods

CAD model of
hip implant




ASTM F2182
phantom




Birdcage coil
resonant at
64MHz (1.5T)

Materials and methods



Two implant models



Various parts/materials



FDTD technique



Software used


- Coupled electrothermal problem
- Results normalized to 2W/kg in phantom (normal mode)

Temperature increase of tissue due to hip implants inside an MRI birdcage coil


V. Gkyrinis and T. Samaras, Aristotle University of Thessaloniki (GR)

Results

Electric field distribution



Induced electric currents




Temperature increase of tissue due to hip implants inside an MRI birdcage coil


V. Gkyrinis and T. Samaras, Aristotle University of Thessaloniki (GR)

Results

Peak spatial SAR_{1g}



Temperature rise




Temperature increase of tissue due to hip implants inside an MRI birdcage coil

V. Gkyrinis and T. Samaras, Aristotle University of Thessaloniki (GR)

Results

Temperature rise at various points with different materials



Top Implant								
Point :	1	2	3	4	5	6	7	8
Co-Cr-Mo	11.3	3.4	0.9	0.6	0.7	0.9	0.8	1.2
Stainless Steel 316L	11.0	3.3	0.9	0.6	0.7	0.9	0.8	1.1
Ti Alloy Grade 5	12.7	3.5	0.7	0.5	0.7	1.0	0.9	1.5
Steel + Alumina Liner	11.0	3.3	0.9	0.7	1.2	1.7	1.5	1.6
Steel + Zirconia Liner	11.1	3.3	0.9	0.7	1.0	1.5	1.3	1.7
Steel + Polymer Liner	10.9	3.2	0.9	0.8	1.2	1.9	1.6	2.1
Bottom Implant								
Point :	1	2	3	4	5	6	7	8
Co-Cr-Mo	8.2	2.5	1.7	1.6	1.7	3.3	2.6	3.2
Stainless Steel 316L	8.0	2.5	1.7	1.6	1.7	3.2	2.5	3.1
Ti Alloy Grade 5	9.0	2.5	1.5	1.7	1.6	3.5	2.6	3.5
Steel + Alumina Liner	8.3	2.6	1.8	1.8	2.3	4.2	3.6	3.8
Steel + Zirconia Liner	8.3	2.6	1.8	1.7	2.2	4.1	3.3	3.9
Steel + Polymer Liner	8.2	2.5	1.8	1.8	2.4	4.5	3.7	4.2

Conclusions

- Modern numerical tools allow for the prediction of temperature rise in tissue next to an implant inside the MRI
- The comparison of measurements in the literature with present numerical results shows very good agreement
- Computer simulations can be used for parametric analysis of tissue heating with implant shape, material used, etc.
- Costly, lengthy and tedious measurements in an experimental setup could be replaced in the future by calculations